


② single use
sterile packed



nxt™ single-use non-stick bipolar forceps



nxt™ non stick single-use bipolar forceps

Superior non-stick technology - 100% made in Germany

- Solid silver tips for a superior non-stick effect (compared to silver coated tips)
- Slim forceps profile, thin tines and excellent visualization
- Can be used with standard bipolar cables
- Optimal balance and weight
- High stability (no scissoring)
- Precision coagulation
- True tip alignment
- Premium look & feel


nxt™ Non-Stick Single-Use Bipolarpinzetten

Überragende Non-Stick Eigenschaften - 100% made in Germany

- Solide Silberspitzen für einen überlegenen Non-Stick Effekt (im Vergleich zu beschichteten Pinzetten)
- Schmales Pinzettenprofil, dünne Schenkel, exzellente Sicht
- Verwendbar mit Standard Bipolarkabeln
- Optimale Balance und Gewicht
- Hohe Stabilität (kein Kreuzen der Pinzettenspitzen)
- Präzise Koagulation
- Exakt ausgerichtete Spitzen
- Exzellente Verarbeitung und Produktqualität



nxt™
Non-Stick Technology VS. **standard** bipolar forceps



Special nxt™ material

Active heat transport

No sticking of the forceps!

Standard material

Forceps tips overheat

Forceps begin to stick



Rupture of coagulated tissue



wide handle



single-use classic non-stick



Ø 0.4 x 6 mm micro tip straight

standard series / working / total length 50 / 160 mm

BIP047USSU us pin connector

alternative to Spetzler® Malis™ or VersaTru™



Ø 0.7 x 6 mm micro tip straight

standard series / working / total length 50 / 160 mm

BIP048USSU us pin connector

alternative to Spetzler® Malis™ or VersaTru™



Ø 1.5 x 8 mm standard tip straight

standard series / working / total length 50 / 160 mm


BIP048-1.5USSU us pin connector

alternative to Spetzler® Malis™ or VersaTru™

wide handle



single-use classic non-stick




Ø 0.4 x 6 mm micro tip straight

standard series / working / total length 70 / 175 mm

BIP050USSU us pin connector

alternative to Spetzler® Malis™ 6760-180-005 or VersaTru™ 9007050ST




Ø 0.7 x 6 mm micro tip straight

standard series / working / total length 70 / 175 mm

BIP051USSU us pin connector

alternative to Spetzler® Malis™ or VersaTru™




Ø 1.0 x 8 mm standard tip straight

standard series / working / total length 70 / 175 mm

BIP052USSU us pin connector

alternative to Spetzler® Malis™ 6760-180-010 or VersaTru™ 9007100ST




Ø 1.5 x 8 mm standard tip straight

standard series / working / total length 70 / 175 mm

BIP052-1.5USSU us pin connector

alternative to Spetzler® Malis™ 6760-180-015 or VersaTru™ 9007150ST

single-use classic non-stick




Ø 0.2 x 6 mm micro tip straight

standard series / working / total length 95 / 205 mm

BIP053USSU us pin connector

alternative to Spetzler® Malis™ or VersaTru™




Ø 0.4 x 6 mm micro tip straight

standard series / working / total length 95 / 205 mm

BIP054USSU us pin connector

alternative to Spetzler® Malis™ 6760(80)-200-005 or VersaTru™ 9008050ST(SL)




Ø 1.0 x 8 mm standard tip straight

standard series / working / total length 95 / 205 mm

BIP056USSU us pin connector

alternative to Spetzler® Malis™ 6760(80)-200-010 or VersaTru™ 9008100ST(SL)




Ø 1.2 x 8 mm standard tip straight

standard series / working / total length 95 / 205 mm

BIP057USSU us pin connector

alternative to Spetzler® Malis™ or VersaTru™




Ø 1.5 x 8 mm standard tip straight

standard series / working / total length 95 / 205 mm

BIP056-1.5USSU us pin connector

alternative to Spetzler® Malis™ 6760(80)-200-015 or VersaTru™ 9008150ST(SL)




Ø 2.0 x 8 mm standard tip straight

standard series / working / total length 95 / 205 mm

BIP058USSU us pin connector

alternative to Spetzler® Malis™ or VersaTru™




Ø 1.2 x 8 mm standard tip 45° down

standard series / working / total length 95 / 205 mm

BIP059USSU us pin connector

alternative to Spetzler® Malis™ or VersaTru™



Ø 1.2 x 8 mm standard tip 45° up

standard series / working / total length 95 / 205 mm

BIP060USSU us pin connector

alternative to Spetzler® Malis™ or VersaTru™

wide handle




1:1

wide handle



single-use classic non-stick




Ø 0.4 x 6 mm micro tip straight

standard series / working / total length 130 / 235 mm

BIP062USSU us pin connector

alternative to Spetzler® Malis™ 6760(80)-230-005 or VersaTru™ 9009050ST(SL)




Ø 1.0 x 8 mm standard tip straight

standard series / working / total length 130 / 235 mm

BIP064USSU us pin connector

alternative to Spetzler® Malis™ 6760(80)-230-010 or VersaTru™ 9009100ST(SL)




Ø 1.2 x 8 mm standard tip straight

standard series / working / total length 130 / 235 mm

BIP065USSU us pin connector

alternative to Spetzler® Malis™ or VersaTru™



Ø 1.5 x 8 mm standard tip straight

standard series / working / total length 130 / 235 mm


BIP064-1.5USSU us pin connector

alternative to Spetzler® Malis™ 6760(80)-230-015 or VersaTru™ 9009150ST(SL)


wide handle



single-use irrigation non-stick



Ø 0.5 x 6 mm micro tip straight
standard series / working / total length 95 / 205 mm
BIP002US-IRSU us pin connector
alternative to Spetzler® Malis™ 6770-200-005



Ø 1.0 x 8 mm standard tip straight
standard series / working / total length 95 / 205 mm
BIP003US-IRSU us pin connector
alternative to Spetzler® Malis™ 6770-200-010

wide handle



single-use irrigation non-stick



Ø 0.5 x 6 mm micro tip straight
standard series / working / total length 130 / 235 mm
BIP004US-IRSU us pin connector
 alternative to Spetzler® Malis™ 6770-230-005



Ø 1.0 x 8 mm standard tip straight
standard series / working / total length 130 / 235 mm
BIP005US-IRSU us pin connector
 alternative to Spetzler® Malis™ 6770-230-010



adeor® is a registered trademark of adeor medical AG.

All rights of distribution, also by photocopy, reprint or storing and recovery in any kind of data processing device are reserved and require our written approval.

Copyright 2021 adeor medical AG

01-1121-D/E

Printed in Germany

November 2021

Spetzler and Malis are trademarks of Stryker Inc.

VersaTRU is a trademark of Johnson & Johnson Inc.

adeor is not affiliated with any of the above.



HEADQUARTERS

adeor medical AG

Martinshof 5

D-83626 Valley

Germany

Phone +49 (0) 8024 47717 0

Fax +49 (0) 8024 47717 99

eMail office@adeor.com

web www.adeor.com

Printed on recycled paper.

Gedruckt auf Recyclingpapier.

