



ESSENTIAL

**PRESSURECARE
PRODUCTS**





Trulife Overview

Trulife is an internationally based and managed group engaged in the creation, development, manufacture and marketing of niche healthcare products.

In addition to Pressurecare pads for the Operating Table, the Group activities encompass product ranges, Orthopaedics, Prosthetics, Breastcare and Wheelchair cushions. For further details on all our product ranges please visit our website: **www.trulife.com**

The Pressurecare Story

Trulife began making Breast Prostheses in 1987 to make life for women following breast cancer, a little easier. We enjoy a well-deserved reputation as manufacturing the most natural products on the market today.

The knowledge that silicone is a well-recognised material in the medical industry, compatible with skin, durable, and excellent at relieving and re-distributing pressure, encouraged the Trulife team to pursue other uses for silicone.

In 1995 Trulife began manufacturing pressure relieving products for the operating theatre, helping prevent patients from developing pressure sores and nerve injury during surgery. We are now market leaders in many countries around the world.

Mission Statement

The core philosophy is a belief in serving customers' needs through excellence in product innovation, quality in everything we do, extensive education and training offerings, and development of patient care through committed and dedicated people.



Oasis Elite Closed Head Ring

Unique gel & foam combination conforms comfortably to the patient's head, this head ring is suitable for use in many surgical procedures.



Code	Description	Dimensions
EL012	Adult	200 x 50mm



Oasis Elite Prone Head Rest


Developed in conjunction with the Operating Staff at a leading hospital, the Oasis Elite Prone Head Rest is anatomically designed so that the patient's head can comfortably rest face down with respiratory tubes in place. The unique design provides excellent support, stability, comfort and protection during surgery in the prone position, with openings on both sides. The product is lightweight, durable and easy to handle.

Code	Description	Dimensions
EL141	Adult	280 x 240 x 140mm



Oasis Elite Open Head Ring

Combination of gel and foam provides stability and protection to the patient's face during surgery. An opening allows for intubation.



Code	Description	Dimensions
EL022	Adult	200 x 54mm



Oasis Elite Supine Head Rest

This Supine Head Rest is ideal for positioning the patient during eye surgery and other supine procedures. Ergonomically shaped to position the patient's head, while effectively protecting against the onset of pressure sores.

Code	Description	Dimensions
EL151	Adult	275 x 237 x 90mm



Oasis Elite Contoured ArmBoard Pad

This anatomically shaped armboard pad is designed to provide added support and stability to the arm in a variety of procedures.



Code	Description	Dimensions
EL401	Adult	510 x 155 x 40mm



Oasis Elite Heel Pads

The Oasis Elite Heel Pads provide support and maximum protection for the heels, even during the longest surgical procedures. They support the Achilles and effectively remove pressure from the heels. Supplied in pairs, the products are lightweight and easy to use.



Code	Description	Dimensions
EL130	Adolescent/Adult	180 x 100 x 70mm



Oasis Elite Flat Bottom Chest Rolls

Chest Rolls provide excellent support and stability. The unique gel and foam combination means they are up to 75% lighter than comparative gel products, allowing for easy manipulation and handling by OR staff.

Code	Description	Dimensions
EL216	Adult	505 x 140 x 110mm
EL218	Adult	510 x 150 x 150mm




Oasis Crutch/Universal Pads

The Oasis Crutch Stirrup Pads are designed to fit most types of stirrups. Supplied in pairs, they come in a variety of sizes.



Code	Description	Dimensions
OA091	Adolescent/Adult	395 x 245 x 10mm


Oasis Universal Patient Positioners



These 100% silicone positioners are suitable for safe and effective positioning during all surgical procedures.

Code	Description	Dimensions
OA053	Adolescent/Adult	345 x 120 x 70mm

Oasis Standard Operating Table Pads



These Oasis table pads provide protection over the entire surface area of the operating table and are recommended for use with patients in all positions. The pads are flexible, soft, comfortable and versatile and distribute weight evenly over the entire surface area.



Code	Description	Dimensions
OA030	Universal Table Pad	520 x 520 x 10mm
OA031	3/4 Length Table Pad	1150 x 520 x 10mm
OA032	Full Length Table Pad	1800 x 520 x 10mm



Trulife Group Headquarters

Airton Road, Tallaght, Dublin D24 R859

Tel: +353 (0)1 451 1755

www.trulife.com